


## Product specifications

### SpinChem® Vessel V2

2221-001

#### Properties

<b>Material:</b>	Borosilicate glass
<b>Height:</b>	250 mm
<b>Width:</b>	140 mm
<b>Weight:</b>	763 g
<b>Nominal volume:</b>	200 mL
<b>Flange:</b>	DN100
<b>Hose connection:</b>	GL18
<b>Lid height:</b>	110 mm (155 mm including shaft guide)
<b>Lid weight:</b>	549 g
<b>Lid necks:</b>	5 ea; straight B24, angled B34, B29, B24, B19



#### Operating conditions

<b>Pressure:</b>	0 to 0.5 bar
<b>Temperature:</b>	-20 to +110 °C*
	*-70 to +230 °C without hose connectors
<b>Liquid volume:</b>	120 to 300 mL*
	*Determined with water at 20°C and a fully packed RBR S2 rotating at 500 rpm