

Product specifications

SpinChem® RBR S2

1221-001

Properties

Material:	Stainless steel (SS316L / EN2348)
Diameter:	45 mm
Height:	30 mm
Weight:	317 g (85 g without shaft)
Shaft length:	380 mm
Shaft diameter:	10 mm
Volume:	28 mL
Filter porosity:	104 µm



Operating conditions

Rotational direction: Clockwise

Rotational speed: 300 - 1000 rpm

NB: Maximum 500 rpm continuously when shaft guide is used!