

Product specifications

SpinChem® RBR S3

1321-001

Properties

Material:	Stainless steel (SS316L / EN2348)
Diameter:	70 mm
Height:	30 mm
Weight:	397 g (148 g without shaft)
Shaft length:	410 mm
Shaft diameter:	10 mm
Volume:	69 mL
Filter porosity:	104 µm



Operating conditions

Rotational direction: Clockwise

Rotational speed: 200 - 800 rpm

NB: Maximum 500 rpm continuously when shaft guide is used!