

Product specifications

SpinChem® RBR S3 plus

1321-201

Properties

Material:	Stainless steel (SS316L / EN2348)
Diameter:	70 mm
Height:	63 mm
Weight:	537 g (289 g without shaft)
Shaft length:	410 mm
Shaft diameter:	10 mm
Volume:	138 mL
Filter porosity:	104 µm



Operating conditions

Rotational direction: Clockwise

Rotational speed: 200 to 800 rpm

NB: Maximum 500 rpm continuously when shaft guide is used!