

# Product specification sheet



1411-001

## SpinChem® RBR S4

### Properties

<b>Material:</b>	Stainless steel (SS316L / EN2348)
<b>Diameter:</b>	140 mm
<b>Height:</b>	140 mm
<b>Weight:</b>	2.4 kg without shaft
<b>Shaft length:</b>	Customized
<b>Shaft diameter:</b>	30 mm
<b>Bed volume:</b>	1.1 L
<b>Filter porosity:</b>	100 µm

### Operating conditions

<b>Rotational direction:</b>	Clockwise / Counterclockwise
<b>Rotational speed:</b>	150 to 400 rpm



SpinChem AB

✉ Tvistevägen 48C, 90736 UMEÅ, Sweden

☎ +46 90192501