

Product specification sheet



2221-001

SpinChem® Vessel V2


Properties

Material:	Borosilicate glass
Height:	250 mm
Width:	140 mm
Weight:	763 g
Nominal volume:	200 mL
Flange:	DN100
Glass thread:	GL18
Lid height:	110 mm (155 mm including shaft guide)
Lid weight:	549 g
Lid necks:	5 total: straight B24, angled B34, B29, B24, B19

Operating conditions

Overpressure:	0-0.5 bar above atmospheric
Vacuum:	10 mbar absolute
Temperature:	-70 to +230 °C
Liquid volume:	120 to 300 mL*

* Determined with water at 20°C and a fully packed RBR S2 rotating at 500 rpm



SpinChem AB

✉ Tvistevägen 48C, 90736 UMEÅ, Sweden

☎ +46 90192501